



## *PICC: How to Flush/Keep Patent*

Supplies: 10ml Normal Saline flush  
Alcohol prep pad

Procedure:

- a. Wash hands for 20 seconds
- b. Check the expiration date on the Normal Saline flush
- c. Remove all the air from the Normal Saline syringe
- d. Clean the needle less injection cap with an alcohol pad
- e. Insert the syringe into the end of the injection cap
- f. Unclamp the catheter
- g. Flush the catheter with push/pause technique
- h. Remove the syringe
- j. Clamp the catheter

### *Catheter Information*

My catheter brand is: \_\_\_\_\_

My catheter make is: \_\_\_\_\_

My catheter Lot number is: \_\_\_\_\_

I have a Single/ Dual /Tripple Lumen Catheter (circle one)